


# Kidney Clamps Nierenklemmen

---

## ATRAUMATA



# Kidney Stone Forceps Nierensteinzangen

---



08-189




08-190



08-191



1/4  
08-192



08-193



08-194



1/4  
08-195



1/2


08-189 - 192  
RANDALL  
225 mm




1/2

08-193 - 195  
RANDALL  
190 mm


# Bladder Retractors Blasenwundspreizer



**09-101**  
JUDD-MASSON  
Bladder and suprapubic  
retractor




**09-102**  
MILLIN  
Center bladder for  
retro pubic prostatectomy



**09-103**  
MILLIN  
Bladder retractor, complete  
105X45 mm




**09-104**  
LEGUEU  
Bladder retractor  
Center blade  
80 x 31 mm



**09-105**  
KOCHER  
Bladder spatula  
235 mm


# Trocars Trocare



Luer Lock



**09-106**  
KLIMA ROSEGGER  
35 X 1,8 mm  
Sternal puncture needle



**09-107**  
Universal trocar  
with 4 tips




NELSON


**09-108 - 112**  
Inside diam. Outside diam  
**09-108** 2,2 mm 3,0 mm  
**09-109** 3,0 mm 4,0 mm  
**09-110** 4,0 mm 5,0 mm  
**09-111** 5,0 mm 6,0 mm  
**09-112** 5,5 mm 6,5 mm  
for hydrocele and abdominal puncture



**09-113 - 116**  
Inside diam. Outside diam  
**09-113** 2,2 mm 3,0 mm  
**09-114** 3,0 mm 4,0 mm  
**09-115** 4,0 mm 5,0 mm  
**09-116** 5,0 mm 6,0 mm  
for abdominal puncture,  
with stop-cock and lateral  
drain



**09-117 - 118**  
Inside diam. Outside diam  
**09-117** 6,5 mm 8,0 mm  
**09-118** 9,0 mm 10,0 mm  
Epyema trocar  
BUELEU



**09-119 - 122**  
Inside diam. Outside diam  
**09-119** 3,0 mm 4,0 mm  
**09-120** 3,6 mm 4,6 mm  
**09-121** 4,3 mm 5,3 mm  
**09-122** 5,0 mm 6,0 mm  
Gall bladder trocars  
OCHSNER



Record Cone1  
**09-123**  
PIERCE  
Antrum trocar  
Ø 4,0 mm



Luer - Lock  
Cone  
**09-124**  
COAKLEY  
Antrum trocar  
Ø 3,0 mm